

Tracoe Purofoam – REF 958 & 959 Dressings

Product Description

There are many essential products that play an important role in tracheostomy care. Our assortment of essentials includes neckbands, dressings, products for care & cleaning as well as clothing & protectors and accessories. They fulfill a wide range of user needs and preferences, are easy to use, comfortable, discreet, offer security, and are functional.



| General Information | | Material Specifications | |
|-------------------------------|-------------------------------------------------------------------------------|---------------------------------|------------------------------|
| Product line | Dressings | Material Compress | Polyurethane foam |
| Legal manufacturer | Tracoe Medical GmbH Reichelsheimer Str. 1/3 55268 Nieder-Olm Germany | Material Coating | Polyurethane film / net film |
| Country of origin | Germany | Biocompatibility | EN ISO 10993-1:2009 |
| Handling | | | |
| Shelf life | 3 years | Storage | |
| Maximum Useful Life | Up to 3 days | | |
| Packaging | | Certification | |
| Packaging unit | 10 pcs | Product classification acc. MDR | IIa |
| Dimensions (W x H x L) | *190 x 160 x 80 mm **135 x 135 x 80 mm | CE Marking | |
| Material Primary Packaging | Sterile pouch made of Aluminium - PET | GMDN | 46854 |
| Material Secondary Packaging: | Cardboard | Quality system certified to | EN ISO 13485 |
| | | FDA registration | Exempt/E652368 |

| Product Specifications | | | | |
|--------------------------------|-----------|---------------------------|------------|------|
| Product Name | REF | Dimensions (L x W x H) | Weight (g) | Slit |
| Tracoe Purofoam Compress Large | REF 958* | 9.0 cm x 10.0 cm x 0.5 cm | 12,6 | Yes |
| Tracoe Purofoam Compress Small | REF 959** | 6.5 cm x 6.5 cm x 0.5 cm | 6,0 | Yes |



Tracoe Medical GmbH,
Reichelsheimer Str. 1 / 3,
55268 Nieder-Olm
Germany



Breathing-Speaking-Living

Atos Medical AB
Hyllie Boulevard 17
215 32 Malmö
Sweden
www.atosmedical.com