

Tracoe Vario XL – REF 461

Tracheostomy tube spiral-reinforced, with adjustable neck flange, low pressure cuff, scale and 15 mm connector

Product Description



Tracoe Vario is a soft, flexible tube with an adjustable neck flange. By means of a practical push-button mechanism, the flange can be moved and individually adjusted for each patient. In this way, the correct length can be achieved. A scale ensures the exact positioning of the neck flange in preparation for a tube change. The flange is then locked in place by moving the orange lever upward. Locking with a twist lock is no longer necessary. Tracoe Vario tubes are made of DEHP-free, medical grade PVC and are available with and without spiral-reinforcement. The assortment includes both standard and extra-long versions.



Scope of Supply

- Obturator
- Adjustable neck strap
- Product info card with two detachable labels

General Information

Product line	Vario XL
Legal manufacturer	Tracoe Medical GmbH Reichelsheimer Str. 1/3 55268 Nieder-Olm, Germany
Country of origin	Germany
Colour Coding	REF 461 XL
Handling	  
Shelf life	5 years, from date of production
Maximum Useful Life	29 days, from first use

Packaging

Packaging unit	1 piece
Weight	See Table
Dimensions (W x H x L)	210 x 113 x 50 mm
Material outer packaging:	Cardboard box
Material sterile packaging:	Sterile bag with medical paper and transparent heatsealable PET/PP foil

Material Specifications


Material Outer cannula	Polyvinylchloride
Material Neck Flange	Polysulfone/ TPE
Material Cuff	Polyvinylchloride
Biocompatibility	EN ISO 10993-1



Storage



Certification


Product classification acc. MDR	IIb
CE Marking	
GMDN	35404
Quality system certified to	EN ISO 13485
FDA registration	Clearance Number K203362

Product Range

Size	Order Number	ID	OD	TL (TL Min-TL Max)	OB	BA	CD	Gross weight	Global Trade Item Number
		mm	mm	mm	mm	°	mm	g	
07	REF 461-07	7.0	9.7	64.5 – 97.5	110	100	23.0	206	04035324018970
08	REF 461-08	8.0	11.2	73.5 – 113.5	128	100	28.0	211	04035324018987
09	REF 461-09	9.0	12.3	78.0 – 123.0	134	100	30.0	213	04035324018994
10	REF 461-10	10.0	13.7	84.5 – 129.5	146	100	32.0	217	04035324019007

Legend

ID	Inner diameter at bottom of cannula
OD	Outer diameter at bottom of cannula
TL MIN	Minimum length along centre line
TL MAX	Maximum length along centre line
OB	Length along outer bend
BA	Bending angle
CD	Cuff diameter



Tracoe Medical GmbH
Reichelsheimer Str. 1 / 3
55268 Nieder-Olm
Germany



Breathing-Speaking-Living

Atos Medical AB
Hyllie Boulevard 17
215 32 Malmö
Sweden
www.atosmedical.com